



Solutions **with care**

Pressure Relief & Patient Positioning



SUPINE & LITHOTOMY POSITION

Head Base Ring



- The Head ring gel pads can be placed under the patients head during supine positions
- The unique shape conforms to the patients head and helps in the healing of a sore occipital bone
- The Head Base Ring is available in adult, pediatric and neonatal sizes

Product Code 060A001	Dim (mm) Φ 200 X 45
-------------------------	------------------------

Product Code 025A001	Dim (mm) Φ 140 X 35
-------------------------	------------------------

Product Code 061A001	Dim (mm) Φ 90 X 20
-------------------------	-----------------------

Head Base Gel Pad

- The Head base gel pad can be placed under the patients head during supine positions
- The unique shape conforms to the patients head
- The Head base gel pad is available in neonatal and adult sizes

Product Code 020A001	Dim (mm) Φ 200 X 50
-------------------------	------------------------

Product Code 021A001	Dim (mm) Φ 90 X 20
-------------------------	-----------------------



Head Neck Base Gel Pad



- The Head neck base can be placed under the neck or inferior to the occipital bone of the head to provide extension of the neck muscles during surgeries in the region
- The product can also be used to support other bony prominence especially in the lower hip

Product Code 022A001	Dim (mm) Φ 160 X 38
-------------------------	------------------------

Universal Pad



The Universal gel pad can be used as:-

- Pediatric gel pad for pediatric patients during surgeries
- Incubator pad providing comforting base for neonates from protection against nerve damage
- Hip pad for orthopedic table tops
- Mattress overlay in hospital beds

Product Code
040A001

Dim (mm)
⌀ 510X350X10

Product Code
041A001

Dim (mm)
⌀ 510X350X25

Ulnar Protector

- The Ulnar protector serves to protect the nerves of the elbow and the arm
- The product wraps around the arm and is secured in place with velcro straps
- The product can be used when the patient is positioned in supine, lithotomy or lateral

Product Code
069A001

Dim (mm)
⌀ 460X210X10



Sacral Pad




- The Sacral pad is commonly used for supine and lithotomy positions when there is not enough support to the sacral area
- The Sacral pad provides the required support to the patients not only on the sacral area but also in the peripheral region

Product Code
043A001

Dim (mm)
⌀ 400X380X25

Crutch Gel Pad



- The Crutch gel pads are designed for use during a gynecological procedure, especially when the patient is placed in a lithotomy position
- The Crutch gel pads are designed to fit most of the stirrups in the operating or examination room

Product Code
030A001

Dim (mm)
Φ 320X220X10

Heel Pad


- The Heel cups are used for longer duration surgeries to be placed under the heels of the patient. It protects the heel and the achilles
- The product is provided in pair to prevent any risk of skin breakdown or damage
- The Heel cups are primarily used during the supine position of the patient



Product Code
089A001

Dim (mm)
Φ 190X120X70

Supine Head Rest



- Supine head rest is ideal for positioning the patients head for eye procedures or related supine procedures
- The contour on the supine head rest cradles the patients head during the procedure
- The product is lightweight and provides excellent stability

Product Code
068A001

Dim (mm)
Φ 275 X 240 X 90

PRONE POSITION


Horseshoe



- The Horseshoe gel pad can be used when the patient is in prone position
- The Horseshoe shape provides space for the intubation tubes
- The Horseshoe gel pad is available in adolescent and adult sizes
- The Horseshoe configuration is also available without a center base to be used for similar purposes

Product Code 050A001	Dim (mm) Φ 200 X 50
Product Code 055A001	Dim (mm) Φ 140 X 40
Product Code 051A001	Dim (mm) Φ 90 x 20

Horseshoe Without Base



- The Horseshoe gel pad can be used when the patient is in prone position
- The Horseshoe shape provides space for the intubation tubes

Product Code 050A001 -1	Dim (mm) Φ 200 X 50
----------------------------	------------------------

Split Horseshoe



- The Split horseshoe provides cranial support in prone position
- The Split horseshoe can be provided with and without a frame that is compatible with Mayfield and Sugitha headrest
- The Split horseshoe can be attached to the frame by inserting the metal horseshoe frame inside the cloth pocket of the split horseshoe
- The Split horseshoe gel pad is available in adolescent and adult sizes

Product Code 050A001-2	Dim (mm) Φ 180x45x50
Product Code 055A001-2	Dim (mm) Φ 115x35x40

Split Horseshoe With Frame



- The frame allows for lateral and vertical adjustments and hence flexibility in patient positioning
- The width of the horse shoe can be adjusted to accommodate various head sizes of the patient
- The water resistant cloth of the gel Split horseshoe enables easy cleaning of the product after use
- The Split horseshoe gel pad is available in adolescent and adult sizes

Product Code 050A001-3	Dim (mm) Φ 180x45x50
Product Code 055A001-3	Dim (mm) Φ 115x35x40

Flat Prone Gel Bed



- The Flat prone gel bed provides very good support to the patient placed in a prone position
- It eliminates the need to use the bolsters which often require accurate positioning while the patient is placed in prone
- The Prone gel bed design relieves the abdomen thereby assisting the patient in keeping a normal breathing rate
- The Flat prone gel is commonly used for spinal procedures

Product Code
053A001

Dim (mm)
Φ 650X460X130

X-ray Translucent
Ideal for spinal procedures
Facility to adjust width of the bed.

X-ray Translucent
Ideal for spinal procedures
Facility to adjust width of the bed.

The Prone position gel bed can be used to

- Relieve the chest in prone position thereby maintaining better respiratory function
- Protect the hip and chest for extended surgeries requiring prone position

Product Code
052A001

Dim (mm)
Φ 600X490X180

Prone Position Gel Bed



Knee/Arm Rest Gel Pad



The Knee and the Arm-board gel pad can be used to protect

- Knee and the Arm in prone position
- Upper arm in prone and lateral positions
- The Knee/Arm rest gel pad is available in adolescent and adult sizes

Product Code
070A001

Dim (mm)
Φ 500X150X35

Product Code
070A001-1

Dim (mm)
Φ 350X150X35

Prone Head Rest Gel



- The Gel prone head rest is developed to support the patients head during surgeries requiring the patient to be in prone position
- The Gel prone head rest design supports the forehead and the chin of the patient
- The Prone head rest gel does not cause any pressure on the eyes of the patient that may result in an increase in the intraocular pressure
- There is sufficient space for the anesthesia tubes to pass through without kinking
- The Gel prone head rest is an ideal choice for longer duration surgeries

Product Code
056A001P

Dim (mm)
Φ 270X220X135

Dome Positioner Prime

The Dome positioner or also the chest roll can be used to:

- Support shoulder during spine surgery
- Support the back during lateral positioning
- Support under the knee when required
- To raise and protect the heels when required

Product Code 080A001P	Dim (mm) Φ 510X150X160
Product Code 081A001P	Dim (mm) Φ 510X150X120
Product Code 087A001P	Dim (mm) Φ 510X150X90
Product Code 087A001- 1P	Dim (mm) Φ 350 X 150 X 90



LATERAL POSITION

Wedge Lateral Positioner




- The Lateral positioner is a Wedge shaped gel positioner that can be used to position the patient in a semi lateral position or semi prone position
- In both the positions, the Wedge shaped lateral positioner can be placed either on the back or near the chest respectively

Product Code
065A001

Dim (mm)
Φ 315X265X100

Lateral Gel Positioning Bed


- The Lateral positioning system allows the anesthetist to position the patient in a lateral position from a supine position safely and effectively
- The U shaped cradle allows the patients shoulder to fit in it hence providing comfort and nerve protection
- The Lateral positioning system is used for shoulder surgery, hip replacement, kidney and other related surgical procedures



Product Code
066A001

Dim (mm)
Φ 720x500x160

Lateral Leg Positioner - Foam with Gel Top



- The Lateral leg positioner is used to position the legs of the patient in the lateral position
- The Lateral leg positioner is available with and without a gel top
- While using the product when the patient is in lateral positioner, the top leg is placed on the gel portion while the other leg is underneath
- The positioner helps to maintain hip alignment and relieve pressure on the knees during procedures

Product Code
067A001

Dim (mm)
Φ 650x400x240

Product Code
067A001-1

Dim (mm)
Φ 650x400x230

OT TABLE GEL MATS

OT Table Full Length Gel Pad



- Can be placed on the operating table for most surgeries
- Prevent nerve damage and better protection than regular foam mattress
- Are available in two standard sizes, full length and pediatric

Product Code 057A001	Dim (mm) Φ 1750x510x10
-------------------------	---------------------------

Product Code 057A001-1	Dim (mm) Φ 1050x510x10
---------------------------	---------------------------

OT Table Gel Pad - Head Section/ Torso Section/ Leg Section

- The three piece gel pad sections, i.e. Head section, Torso section and Leg section can be used for most surgical procedures
- The sectional gel pads offer versatility when using them during surgery
- The sectional gel pads prevent nerve damage and are available in two thicknesses; 10mm and 25 mm
- The thicker gel pads offer higher degree of protection and comfort

Product Code 058A001	Dim (mm) Φ 220x500x25
-------------------------	--------------------------


Product Code 058A001- 1	Dim (mm) Φ 900x500x25
----------------------------	--------------------------

Product Code 058A001- 2	Dim (mm) Φ 500x500x25
----------------------------	--------------------------



COMMON GEL PADS

Flat Gel Pad 15 & 35




- The Flat gel pad is a versatile product that can be used at various locations on the body where the anesthetist or the surgeon feels that pressure relief is required
- The Flat gel pad comes in two different thicknesses
- The pad can also be used under the head, shoulders, and sacrum
- The choice of the thickness will depend on the position and the area that is being protected

Product Code 063A001	Dim (mm) Φ 260X220X15
-------------------------	--------------------------

Product Code 064A001	Dim (mm) Φ 260X220X35
-------------------------	--------------------------

Flat Gel Pad 35 with Center Pocket

- The Flat gel pad 35 with center pocket is a versatile product that can be used at various locations on the body where the anesthetist or the surgeon feels that pressure relief is required
- The pad can also be used under the head, shoulders, sacrum and other possible areas
- The center pocket provides relief to a specific bony prominence that can otherwise not be supported or relieved with other products



Product Code 064A001-1	Dim (mm) Φ 260X220X35
---------------------------	--------------------------

Coccydynia Cushion with Gel Top - Hard/ Soft



- Coccydynia cushion is made of foam or gel with a cutout at the back near the base of the spine to relieve pressure on the tailbone or spine
- These cushions are also made with gel tops. Gel top layer can help provide a softer form of cushioning and will conform better to the contours of your body
- They are used to treat a variety of spinal and back pain disorders
- The cushion can be used on a chair with a back for extra support
- The cushion naturally helps improve posture by slightly elevating the hips hence helping the user to sit up straight and eliminate pressure on spine and pelvis

Product Code 090A001-P	Dim (mm) ⌀ 430X410X85
---------------------------	--------------------------

Product Code 090A001-T	Dim (mm) ⌀ 430X410X75
---------------------------	--------------------------

Product Code 090A001-S	Dim (mm) ⌀ 430X410X75
---------------------------	--------------------------

Gel Cushion



- Gel cushion is available in both as a Gel pad and with a thick foam base
- Allows for easy blood circulation, reduces fatigue associated with long term immobility, and provides comfort to the user when seated for longer duration
- Can be used in a hospital setting, home, office and car
- Is covered with a water resistant cloth to allow easy cleaning and light washing
- Is available in single size; can be used by adults and elder children irrespective of weight and size

Product Code 092A001	Dim (mm) ⌀ 430X350X15
-------------------------	--------------------------

Product Code 091A001	Dim (mm) ⌀ 430X350X40
-------------------------	--------------------------

Universal Gel Pillows



- The Universal gel pillows can be used in various surgical positions to help protect the patient during surgery
- The Gel pillows are available in different sizes to be used for adults and pediatric population
- The Gel pillows can replace other methods of positioning the patient such as bean bags in the operating room

Product Code 044A001	Dim (mm) Φ 160X80X40
Product Code 045A001	Dim (mm) Φ 330X70X50
Product Code 046A001	Dim (mm) Φ 330X120X60
Product Code 047A001	Dim (mm) Φ 360X120X70
Product Code 048A001	Dim (mm) Φ 500X120X70

WC Seat Cover


- The Commode seat cover is available in a standard size
- Protects the sensitive nerves in the upper leg and gluteal areas of the patient
- The product can be washed after every use

Product Code 031A001	Dim (mm) Φ 420X400X10
-------------------------	--------------------------



NEONATAL

Neonatal Roll



- The Neonatal rolls can be used to support the neck in the supine position; support and protect the chest of the neonate in the prone position
- The Neonatal rolls is available in various sizes and can be used accordingly for both preterm and full terms

Product Code 082A001	Dim (mm) Φ 38 X 38 X 120
-------------------------	-----------------------------

Product Code 083A001	Dim (mm) Φ 30 X 20 X 120
-------------------------	-----------------------------

Product Code 084A001	Dim (mm) Φ 25 X 25 X 120
-------------------------	-----------------------------

Pediatric Roll

- The Pediatric roll can be used to support the neck in the supine position; support and protect the chest of the neonate in the prone position
- The Pediatric roll is available in various sizes and can be used accordingly for both preterm and full terms



Product Code 082A001-1	Dim (mm) Φ 50X50X180
---------------------------	-------------------------

Neonatal Head Cradle - Preterm/ Full Term



- The Neonatal head cradle is used to support the neonates head while in the warmer or incubators, especially when the neonate is on ventilator or other respiratory support
- It is available in two sizes; preterm and full term size

Product Code 023A001	Dim (mm) Φ 135X150X52
-------------------------	--------------------------

Product Code 024A001	Dim (mm) Φ 185X155X72
-------------------------	--------------------------

Neonatal Full Size Phototherapy Gel Pad



- The Phototherapy gel pad can be used in a neonatal intensive care unit
- The gel pad is transparent and can be used under the infant during double surface phototherapy
- The gel pad protects the delicate skin and nerves of the infant, especially the pre-term infants
- The Phototherapy gel pad can also be used while the infant is placed in the warmer. The gel pad will heat only up to the warmer set temperature and provide the infant with the warmth during the on-off cycles of a standard warmer as well as during an abrupt power disruption
- The Phototherapy gel pad can be washed and cleaned using standard disinfectants



Product Code
042A001 -2

Dim (mm)
Φ 510X350X10

Product List

SUPINE & LITHOTOMY POSITION		PRONE POSITION	LATERAL POSITION	OT TABLE GEL MATS	COMMON GEL PADS	NEONATAL
Head Base Ring	Horseshoe	Wedge Lateral Positioner	OT Table Full Length Gel Pad	Universal Gel Pillow Mini	Neonatal Roll	
Head Base	Horseshoe without Base			Universal Gel Pillow Small	Neonatal/ Pediatric Roll	
Supine Head Rest	Split Horseshoe	Lateral Gel Positioning Bed	OT Table Gel Pad Head Section	Universal Gel Pillow Regular	Neonatal Head Cradle Preterm	
Head Neck Base	Split Horseshoe with Frame			Universal Gel Pillow Large	Neonatal Head Cradle Full Term	
Ulnar Protector	Prone Head Rest Gel	Lateral Leg Positioner Foam with Gel Top	OT Table Gel Pad Torso Section	Universal Gel Pillow Extra Large	Neonatal Full Size Phototherapy Gel Pad	
Sacral Pad	Prone Position Gel Bed			Flat Gel Pad 15		
Heel Pads	Flat Prone Gel Bed	OT Table Gel Pad Leg Section	Flat Gel Pad 35	Flat Gel Pad 35		
Crutch Pads	Dome Positioner			Flat Gel Pad 35 with Center Pocket		
Universal Pad	Knee/ Arm Rest Gel Pad		WC Seat Cover	WC Seat Cover		
				Coccydynia Cushion with Gel Top		
				Coccydynia Cushion Hard		
				Coccydynia Cushion Soft		
				Gel Cushion with Foam Base		
				Gel Cushion		



Lenvitz Medical Solutions Pvt. Ltd.

A: Survey No. 46-2, Narhe Road, Ruby Industrial Estate, Opp. Bank of Maharashtra,
Pune 411041, INDIA

Tel: +91 2024699179 | Mob: +91 9552510301 | Email: info@lenvitz.com